

WCT-TS Software Version 5.10 Release Notes (2/11/2021)

What's New

The newest version of our temperature stage software integrates automatic control of the optional temperature control box and light soaking features. This version also added an option allowing users to maintain stage temperature for a set dwell time. These additions expanded the possible options for measurement sequences and can now be customized in an associated measurement sequence creating software to programmatically change the measurement conditions of the tester.

Bug Fixes and UI Improvements

- Users can now measure at a single point elevated temperature
- Easy result file scrolling
- TCP automation functionality for temperature scans